


## Quantitative Cluster Analysis Applications for N-STORM

Optical super-resolution techniques such as Stochastic Optical Reconstruction Microscopy<sup>1</sup> (STORM) enable the visualization of sub-cellular structures below the diffraction limit, and with high molecular specificity. In this application note we highlight the work of Ricci *et al.*<sup>2</sup>, who have applied STORM microscopy towards probing the distribution of nucleosomes on chromatin fibers, helping to elucidate mechanisms of chromatin folding, gene expression, and pluripotency.

### Using N-STORM to quantify the nano-scale distribution of nucleosomes in intact cell nuclei

The classic view of nucleosome organization is relatively simple: a regular pattern of nucleosomes and histone H1 along a stretch of DNA and forming a fiber about 10 nm in diameter and following the ‘beads on a string’ model. The 10 nm fiber is thought to fold into a fiber about 30 nm in diameter, and to further condense into increasingly higher-ordered structures. However, Ricci *et al.*<sup>2</sup> have applied both single and multi-color STORM imaging to show that nucleosomes organize into heterogeneous nanodomains (Fig. 1c), termed ‘clutches’ in reference to the variably sized egg clutches of birds. Clutch size, composition, and organization are highly variable within and between populations of different cell types.

In this work Nikon’s N-STORM system is demonstrated with ~20 nm lateral resolution, an order of magnitude higher than dictated by the diffraction limit, and thus theoretically capable of resolving 30 nm wide chromatin fibers. However, such fibers have been difficult to identify using conventional optical and electron microscopy techniques. Ricci *et al.*<sup>2</sup> calls the strict existence of the 30 nm fiber into question. Conventional optical microscopies cannot resolve such nano-scale features, while electron microscopies



John R. Allen

Applications and Marketing Specialist, Nikon Instruments Inc.  
1300 Walt Whitman Road, Melville NY 11747-3064, USA

lack the necessary molecular specificity and require harsh sample preparation techniques. Thus STORM is able to fill a unique imaging niche, bridging the gap between resolution and molecular specificity. The work of Ricci *et al.*<sup>2</sup> marks the first time the organization of nucleosomes in intact nuclei has been resolved.

### Using STORM to quantify nucleosome distribution in pluripotent cells


Pluripotent mouse embryonic stem cells (mESCs) were chosen as a model system for evaluating nucleosome clutch size, density, and distribution under several different experimental conditions:

- Treatment with cytokine leukemia inhibitory factor (sLif), a differentiation inhibitor (Fig 2a, b).
- Treatment with inhibitors of Mek and Gsk3 kinases (2iLif), which keeps the cells in a ground state of pluripotency (Fig. 2c).
- A Tcf3 double knockout line (mESC<sup>Tcf3<sup>-/-</sup></sup>). Tcf3 is an effector of the Wnt/ $\beta$ -Catenin pathway, which helps maintain ground state pluripotency (Fig. 2d).
- mESCs differentiated into neural progenitor cells (NPCs), (Fig. 2e).
- An H1 'triple knockout': mESCs with deletion of three H1 isoforms (mESC<sup>H1<sup>tko</sup></sup>), (Fig. 2f).


As expected, the number of localizations per nanodomain and their nnds are lower in pluripotent mESCs with respect to mNPCs. Additionally, Ricci *et al.*<sup>2</sup> demonstrate high anti-correlation between nucleosome clutch size and pluripotency in human induced pluripotent stem cells (hiPSCs) using established gene card technology to score pluripotency using five markers: Nanog, TRA1-60, SSEA4, Oct4, and Sox2.

### Multicolor STORM imaging of H2B with H1 and RNA Polymerase II in a variety of cell lines

2-channel STORM imaging with Nikon's N-STORM system revealed that RNA Polymerase II preferentially associates with small clutches while histone H1 is enriched in large clutches (Fig. 3). This is consistent



with the expectation that smaller and less compacted clutches are more open to the binding of transcription factors and polymerases such as PolII. H1 is important in nucleosome compaction and is enriched in heterochromatin. Furthermore H1 is found to be most abundant in the nuclear periphery, where heterochromatin is most enriched.



## Conclusions

Nikon's N-STORM platform has proven to be a powerful platform for the visualization of multi-protein complexes, as demonstrated by Ricci *et al.* Nikon's new N-STORM 4.0 includes integrated analysis tools for cluster analysis and distance measurements, as well as improved 3D and speed performance. To learn more about Nikon's super-resolution solutions, please visit our website ([www.nikoninstruments.com/sr](http://www.nikoninstruments.com/sr)).

1. Rust, M.J., Bates, M. & Zhuang, X. Sub-diffraction-limit imaging by stochastic optical reconstruction microscopy (STORM). *Nat. Methods* **3**, 793–795 (2006).
2. Ricci, M.A., Manzo, C., Garcia-Parajo, M.F., Lakadamyali, M. & Cosma, M.P. Chromatin Fibers Are Formed by Heterogeneous Groups of Nucleosomes In Vivo. *Cell* **160**, 1145–1158 (2015).